



PRODUCT

Blue Disposable Nitrile Gloves

TECHNICAL DATASHEET

DESCRIPTION

Blue disposable nitrile gloves are single-use, blue, non-sterile and naturally antistatic for use as a physical barrier, particularly to fluids in areas where the user or anyone they come into contact with may have a latex allergy.

Nitrile has a greater resistance to certain chemicals and a higher tensile strength providing a more hard wearing glove. They contain NO natural rubber latex and are powder free. Their applications include use in electronics, general industry and any areas where users have a latex allergy.

TECHNICAL SPECIFICATIONS

Blue disposable nitrile gloves are manufactured in facilities certified against ISO 9001, and tested in accordance with EN 455 Medical gloves for single use Parts 1, 2 and 3, for absence of leaks, physical properties listed above and extractable protein content.

The product is a CE marked class 1 medical device. Gloves intended for medical use may also be referred to as "AQL 1.5". This means that they conform to the minimum level of acceptable physical defects for examination gloves.

FEATURES

- 100% nitrile disposable glove
- Powder free and contains no natural rubber products, naturally antistatic
- A good all round disposable glove

CHEMICAL COMPOSITION

- Acrylonitrile-butadiene
- Ammonia solution
- Sulphur
- Zinc oxide
- Zinc diethyldithiocarbamate
- Tamol
- Ranox LC
- Phenolic antioxidant
- Titanium dioxide
- Coagulant WS
- Colanyl Blue Pigment

These gloves contain NO natural rubber latex or silicon.

To request a quotation or for more information, please call **+1 512-580-4220**
email sales@antistat.com or visit www.antistat.com

IMPORTANT: This data sheet and its contents (the "Information") belong to Antistat or are licensed to it. No licence is granted for the use of it other than for information purposes in connection with the products to which it relates. No licence of any intellectual property rights is granted. The Information is subject to change without notice and replaces all data sheets previously supplied. The Information supplied is believed to be accurate but Antistat assumes no responsibility for its accuracy or completeness, any error in or omission from it or for any use made of it. Users of this data sheet should check for themselves the Information and the suitability of the products for their purpose and not make any assumptions based on information included or omitted. Liability for loss or damage resulting from any reliance on the Information or use of it (including liability resulting from negligence or where Antistat was aware of the possibility of such loss or damage arising) is excluded. © 2024 Antistat.

STORAGE

Store in a cool, dry place and once open keep away from direct sunlight or fluorescent light. Store boxes separately from powdered gloves and latex gloves.

DISPOSAL

Dispose of safely taking into account hazards to health from the product and from contact or contamination during use.

FEATURES OF NITRILE DISPOSABLE GLOVES

I. Pro - Nitrile gloves contain NO natural rubber latex proteins and as such can be classified “Latex Free”. They give superior protection against penetration and microbial infection when compared to latex gloves. They also provide a secure barrier against a wide range of chemicals and solvents. These features are combined with high durability and an elasticity almost comparable to latex but much better than vinyl, giving reasonable fit and providing the user with sensitivity and dexterity in use.

II. Con - Although the fit and feel of nitrile gloves is greatly improved it is still not as good as latex. Nitrile gloves are, on the whole, more expensive than latex or vinyl alternatives. The absence of latex proteins does not stop the possibility of allergic reactions to the chemicals used in the glove’s manufacture.

SYMBOLS USED

The following symbols are used on the inner and outer packaging.

NB: These are the most up to date symbols, older packaging is being constantly updated.

latex free



powder free



single use



To request a quotation or for more information, please call **+1 512-580-4220**
email sales@antistat.com or visit www.antistat.com

IMPORTANT: This data sheet and its contents (the “Information”) belong to Antistat or are licensed to it. No licence is granted for the use of it other than for information purposes in connection with the products to which it relates. No licence of any intellectual property rights is granted. The Information is subject to change without notice and replaces all data sheets previously supplied. The Information supplied is believed to be accurate but Antistat assumes no responsibility for its accuracy or completeness, any error in or omission from it or for any use made of it. Users of this data sheet should check for themselves the Information and the suitability of the products for their purpose and not make any assumptions based on information included or omitted. Liability for loss or damage resulting from any reliance on the Information or use of it (including liability resulting from negligence or where Antistat was aware of the possibility of such loss or damage arising) is excluded. © 2024 Antistat.

POSSIBLE REACTIONS ARE

a. Delayed hypersensitivity reaction (Type IV hypersensitivity):

Symptoms include: dermatitis / skin itching / blistering - typically restricted to the glove contact area. Onset of symptoms is variable (6 – 48 hrs of exposure) and depends on exposure frequency as well as strength of the sensitising allergen.

b. Irritant Contact Dermatitis:

A non-allergic reaction to residual chemicals, restricted to areas of direct contact with onset symptoms in minutes to hours. This chemical dermatitis can also occur with synthetic gloves, scrubs, soaps, disinfectants and excessive powder. Symptoms disappear once contact has stopped.

GENERAL ADVICE

- Always wet hands with warm water before applying soap
- Always dry hands thoroughly
- Do not wear gloves for long periods of time
- Do change gloves after each activity
- Do not wear gloves unnecessarily

PRODUCT CODE	DESCRIPTION	SIZE	QUANTITY
109-0201	Blue Disposable Nitrile Gloves	Small	Box of 100
109-0202	Blue Disposable Nitrile Gloves	Medium	Box of 100
109-0203	Blue Disposable Nitrile Gloves	Large	Box of 100
109-0204	Blue Disposable Nitrile Gloves	Extra Large	Box of 100

To request a quotation or for more information, please call **+1 512-580-4220**
email sales@antistat.com or visit www.antistat.com

IMPORTANT: This data sheet and its contents (the "Information") belong to Antistat or are licensed to it. No licence is granted for the use of it other than for information purposes in connection with the products to which it relates. No licence of any intellectual property rights is granted. The Information is subject to change without notice and replaces all data sheets previously supplied. The Information supplied is believed to be accurate but Antistat assumes no responsibility for its accuracy or completeness, any error in or omission from it or for any use made of it. Users of this data sheet should check for themselves the Information and the suitability of the products for their purpose and not make any assumptions based on information included or omitted. Liability for loss or damage resulting from any reliance on the Information or use of it (including liability resulting from negligence or where Antistat was aware of the possibility of such loss or damage arising) is excluded. © 2024 Antistat.